

NaI Gamma-Ray Spectroscopy

Thursday, January 11, 2018

9:16 AM

Station 2 Partner: Nick Quartermont

Insert 1-5

4] Pre-Lab Prep

IT = Isomeric Transition

→ Gamma Ray or Internal Conversion

Photo Peak = Full energy deposition

a) Review decay schemas for

^{22}Na

β^+

215.54

keV

835.00

keV

γ

1274

keV

^{22}Na

β^+

$^{22}_{10}\text{Ne}$

γ

^{60}Co

β^-

317.32

γ_1

1173.24

keV

γ_2

1332.508

^{60}Co

β^-

γ_1

γ_2

^{109}Cd

ϵ 100%

γ 100% 88.0336 KeV

X-ray 22.163
21.99

2.634-3.748

24.9118

^{60}Ni
28

^{109}Cd

ϵ

γ

^{109}Ag
47

^{137}Cs

β^- 513.97 94%
1175.63 5%

γ 283.5 ~0%
661.659 94.36%

X 32.1939

^{137}Cs

β^-

X

γ


^{137}Ba
56

b)

Gamma Spectra detectable
NaI (TI) γ and X rays
La Br₃

6] Procedure

LaBr_3



Smaller than the NaI


Both require a High voltage connection

NaI(Tl) 3"

5) Equipment



3x3 NaI(Tl)




Anode

High Voltage

Dynode

2x2
 NaI(Tl)



LaBr_3



Back contains output
for test (hookup to Pulser)



Dac contains output
for test (hookup to Pulser)
and output
→ to Oscilloscope
→ to Linear Amp


Sources Used

	Activity	Time Rec	Source Nom	Uncertainty
^{22}Na				
^{57}Co				
^{60}Co				
^{109}Cd				
^{137}Cs	10.02 μCi 370.7 KBq	1 Aug 98 12 PST	6199-442	3%


6. Procedure

a) Set up 3x3" NaI (Tl)


1. Set up



MCB Hookup
to Laptop




Lie down to prevent
rolling.
Detection from
other side



Pre Amp

HV
to Scintillator





to HV Brass

connect
93Ω terminator



to Preamp

Note for lab
 $100\Omega \approx 93\Omega$
Learning only



2. On Amp

.5 μ s	Shaping time
POS	Input
Uni Polar	Output
Preamp = 200 pF	(capacitance)



increase capacitance

95 Ω Terminator on dynode

Run Gamma Vision

Set Conversion Gain to 1024/
(MCB = 1024 channels/bins)

Acquire


↳ MCB Properties

↳ Conversion Gain

↳ arrow select

Note: 10% change in HV bias can change
the signal gain by a factor of 2

3. Snapped Oscilloscope to read amplifier output



$V_{max} = 1.48 \text{ V}$

Pulsar

3.34

x1 attenuator

Pos Output

Amplifier


FG = .68

CG = 20

Shape = .5 μs

4. Disconnect Pulsar from Pre Amp

Set HV Bias to +1000 V



1 kV

.5

5.00

Note! Make sure everything is at zero when starting


Get 137 Cs Source

TO 89 - used to look up info in Certificate Folder

Note: Trigger very Important causes signal to disappear
Determine wanted signal by seeing it move in accordance to distance

Placed at 25 cm (centered)

4.4 V



Stable the closer it is to the detector

Amplifier

FG = .5


CG = 20

Shape = .5 μ s

Ch 2 = Amplifier

Ch 2 = Pre Amp

25 cm



Amp 4.4 V
Rise 800 ns (90% of Amp)
Fall 1.4 μ s (90% to 10%)

Sumpy, multiple signals

25 cm



Amp 1v

Ignoring since it messed up Amplifier signals
Uses primarily to check for clipping